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Battery charger with 9 V rechargeable battery for HK-pregnancy detectors Operators instruction

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English

This battery charger is a high-quality long-range power pack suitable for AC voltages between 100 and 240 V. This power pack must only be operated in connection with the HK pregnancy testing units for sows, sheep and goats starting with the November 2001 model, and a 9V rechargeable battery. **Please note:** Your HK pregnancy tester must under no circumstances be connected to this battery charger if the non-rechargeable battery has been inserted, since this will result in severe damage.

1. Exchanging the battery

To exchange the battery for the 9V rechargeable battery, proceed as follows: Unscrew the two screws on the battery compartment at the back of the appliance with a Phillips screwdriver. Take off the lid, remove the battery from the clip and insert the rechargeable battery. Close the battery compartment. You must make absolutely sure that the lid is properly fitted. Attention: avoid squeezing the cable!

2. The charging process

Open the blue sealing cap on the back of the appliance (illustration 1) with a coin by turning it 90 degrees to the left (illustration 2) and tilt it sideways. You will see the charging socket (illustration 3). Put the round plug of the battery charger into the socket and connect the mains plug with the power supply (illustration 4). The yellow diode next to the charging socket will illuminate

during the charging process. The HK pregnancy tester must be switched off during charging. Charging takes a maximum of 8 hours, after which the yellow diode extinguishes. If you notice the diode lighting up at a later stage, then this merely indicates the conservation charge which counteracts the running-down of the rechargeable battery. The charging of the rechargeable battery is definitely completed within this time, and the round plug can be unplugged from the charging socket. Close the sealing cap with a 90 degree turn to the right. The appliance is operational again.

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Albert Kerbl GmbH